

No	Medical Device	BRAND	ESTABLISHMENT NAME	MD REGISTRATION NO	Complaint	Test	Reference Standard	Summary of Result	Supported/ Not Supported
1	TRANSFER SET	Peritone Youwei / Lucenxia Transflo	ADVENTA HEALTHCARE SDN. BHD.	GA84078400417	Leakage at twist clamp during use.	Hypothesis 1: Micro-tear at tubing			
						Test 1: Fluid leakage	ECRI Test Method	Leakage at twist clamp due to degradation of silicone tubing from prolong exposure to Povidone Iodine	SUPPORTED
						Hypothesis 2: Insufficient connector bond strength			
						Test 1: Static Tensile 15N for 15s	ISO 80369	No detachment of tubing from connectors	Not Supported
						Hypothesis 3: Poor mechanical properties of tubing			
						Test 1: Tensile	ISO 10555-1	Good elastic properties were observed	Not Supported
						Hypothesis 4: Material composition causing defect			
						Test 1: Fourier Transform Infrared (FTIR)	ECRI Test Method	No tubing material degradation	Not Supported
						Hypothesis 5: Kinking induced tubing leakage			
						Test 1: Kink	ECRI Test Method	No impact on flow rate	Not Supported
Hypothesis 6: Material fatigue from repeated clamping									
Test 1: Repeated fatigue from repeated clamping	ECRI Test Method	No tubing leak	Not Supported						
Additional findings									
Test 1: Visual inspection on MiniCap	ECRI Test Method	Crack on MiniCap	Not Applicable						
Test 2: Visual inspection on tubing	ECRI Test Method	Blockage of flow	Not Applicable						
2	SCALPEL BLADE	Connecx-Surgimecx	ADVENTA HEALTHCARE SDN. BHD.	GB4223221-80686	Scalpel blades become blunt mid-procedure	Hypothesis 1: Poor material quality			
						Test 1: FESEM	ECRI Test Method	No contamination on blade surface	Not Supported
						Test 2: EDX	ECRI Test Method	Presence of typical metallic elements	Not Supported
						Hypothesis 2: Excessive bending or cosmetic defects			
						Test 1: Microscopic inspection	ECRI Test Method	No bending or other cosmetic defects found	Not Supported
		Hypothesis 3: Surface contamination							
		Test 1: Fourier Transform Infrared (FTIR)	ECRI Test Method	No distinct peak from FTIR		Not Supported			
		Additional Finding 1: Simulation of scalpel blade (Incorrect cutting technique & unable to cut test medium)							
		Test 1: Cutting technique & ability to cut a medium	ECRI Test Method	Product testing does not reveal any defects		Not supported			
		Additional Finding 2: Corrosion from exposure to local anesthetics & bodily fluids							
Literature review	N/A	Corrosion highly affects incising capability of scalpel blades & increases the risk of tissue trauma	SUPPORTED						
Additional Finding 3: Lack of coating on blade surface									
Literature review	N/A	Accelerating the blunt process & increases tissue resistance	SUPPORTED						
3	Hypromellose Eye Drop	Ain Medicare's Ophthal Mellose Eye Drop Solution	AIN MEDICARE SDN BHD	GB5836723-114988	Eye irritation, burning, allergic reactions	Hypothesis 1: pH deviation			
						Test 1: pH	USP <791>	Value within the human physiological pH tear range (6.5 to 7.6)	Not supported
						Hypothesis 2: Osmolality			
						Test 1: Osmolality	USP <785>	Value within the human physiological osmolality tear range (270 to 310 mOsm/kg)	Not supported
						Hypothesis 3: Leachable chemicals / chemical contamination			
						Test 1: Fourier Transform Infrared (FTIR)	ECRI Test Method	No contamination	Not supported
						Hypothesis 4: Compromised sterility			
						Test 1: Sterility	ECRI Test Method	No growth	Not supported
						Hypothesis 5: Endotoxin/pyrogen contamination			
						Test 1: Endotoxin	ECRI Test Method	Does not contain significant residual bacterial endotoxins	Not supported
Hypothesis 6: Inappropriate Prescription									
Test 1: Misapplication User Survey Analysis Data	ECRI Test Method	Survey results revealed clear prescribing patterns of Ain Ophthal Mellose Eye Drop to moderate to major eye trauma/injury patients, not adhering to its intended purpose of "for minor irritations"	SUPPORTED						
Hypothesis 7: Combined eye drop application									
Test 1: Analysis of Co- Prescription Medication Survey Data	ECRI Test Method	Eye drop does not react with any of the medicine that are prescribed and used together	SUPPORTED						
Hypothesis 8: Cap Leak									
Test 1: Cap Leak	ECRI Test Method	Demonstrated that the reclosable cap of the Ain Ophthal Mellose ampoule does not maintain an adequate seal after first opening	SUPPORTED						
Hypothesis 9: Compromised Cap Integrity									
Test 1: Re-closable Cap Integrity	ECRI Test Method	Demonstrated defective opening mechanisms (poor tamper point scoring or incomplete severing) in Ain Ophthal Mellose ampoules	SUPPORTED						

No	Medical Device	BRAND	ESTABLISHMENT NAME	MD REGISTRATION NO	Complaint	Test	Reference Standard	Summary of Result	Supported/ Not Supported
4	SURGICAL GLOVES	SURGEA surgical gloves (powdered)	AVECENA GLOVES SDN BHD	GB2884724-165273	Surgical gloves tears easily, contains excessive powder and causes allergies (rashes).	Hypothesis 1: Additive in chemical formulation & degradation under high temperature Test 1: Tensile test Test 2: Tensile strain Hypothesis 2: Elevated protein content levels Test 1: Protein content Hypothesis 3: Elevated powder content levels Test 1: Powder content Hypothesis 4: Donning & doffing Test 1: Donning & doffing of glove Hypothesis 5: Excessive appearance of powder Test 1: Powder appearance after 30 minutes	ASTM D3577 ISO 10282 EN455-3 ASTM D6124 ECRI Test Method ECRI Test Method	Exceeds minimum requirement of tensile strength Exceeds minimum requirement of tensile strain Passes protein content levels of below 50 µg/g Passes powder content levels of 15mg/dm2 No difficulty of donning & doffing Granular powder appearance matches powder content quantitative testing	Not supported Not supported Not supported Not supported Not supported Not supported
5	HEMODIALYSIS CATHETER KIT	REAL FLOW hemodialysis catheter from Haolang Technology	ROCKY SCIENTIFIC INSTRUMENTATION	GB8127224-173483	Flow issues due to kinking of catheter and insertion difficulty	Hypothesis 1: Stiffness of catheter body material Test 1: Tensile Hypothesis 2: Incompatible hospital disinfectant to the material Test 1: Fourier Transform Infrared (FTIR) Hypothesis 3: Soft catheter tips cause bending leading to difficult penetration Test 1: Tensile Hypothesis 4: Inconsistency in catheter's physical parameters Test 1: Dimensional measurement of ID, OD and wall thickness Hypothesis 5: No coating on catheter surface Test 1: Talcum powder Hypothesis 6: Catheter tip design led to difficult penetration Test 1: Dimensional measurement of ID, OD and wall thickness Hypothesis 7: Catheter kinking Test 1: Kink	ISO 10555-1 ECRI Test Method ECRI Test Method ECRI Test Method ECRI Test Method ECRI Test Method ECRI Test Method	Real Flow Young's Modulus comparable with Mahurkar Silicone tubing is compatible with alcohol-based CHLOROHEXIDINE Real Flow catheter tip soft and easy to bend Consistent wall thickness was observed Presence of coating on catheter surface Real Flow catheter tip wall thickness was thinner and shorter taper length compared to alternative brands Real Flow flow reduction about 10% due to catheter kink	Not Supported Not Supported SUPPORTED Not Supported Not Supported SUPPORTED SUPPORTED
6	HEMODIALYSIS FLUID FILTER	BLT Biolight F210 fluid filters (Guang Dong Biolight Meditech Co., Ltd.)	VISION BIOMED SERVICES SDN.BHD	GB4794825-194066	Reduction in flow through heater, infection due to different membrane types enabling bacteria passage and leakage	Hypothesis 1: Membrane pore size too large Test 1: FESEM Hypothesis 2: Poor endotoxin retention capability Test 1: Bacteria & Endotoxin Retention Analysis Hypothesis 3: Fluid filter connection incompatibility Test 1: Functional Test & Inspection Hypothesis 4: Connector dimension mismatch Test 1: Filter design & measurement Additional findings: Technical specification Test 1: Review of Data Sheet	ECRI Test Method ISO 23500-5:2024 ECRI Test Method ECRI Test Method	DF220 has slightly larger pore size and pore area than DIASAFE® Plus but the differences are negligible. Meets ISO 23500-5:2024 requirement. The total bacteria count (TBC) and endotoxin level are below <0.1CFU/ml and 0.03EU/ml. No leakage, machine alarms/ errors observed Small differences in connectors, O-rings, and membrane thickness, with no sealing or leaking issues Noticeable difference in Transmembrane Pressure (TMP)	Not Supported Not Supported Not Supported Not Supported SUPPORTED